

## INTRODUCTION

Bovine viral diarrhoea (BVD) is one of the most important infectious cattle diseases worldwide. In Germany, BVD is a notifiable disease, and a nationwide mandatory eradication program has been in place since 2011. Ear notches are prime samples because their collection is easier thanks to identification of animals with ear tags and simultaneous sampling just after birth. The objective is to remove permanently infected animals (PI animals) from the herd at an early stage. Optimized lab workflows are needed, as laboratories handle high number of samples. Test performance, robustness and accuracy is crucial to avoid missing PI animals. In this context, IDvet developed a new real-time RT-PCR kit, the ID Gene™ BVD / BD Triplex 2.0 kit and a reliable direct lysis protocol which significantly simplifies the workflow in the laboratory while increasing the result accuracy (ID GENE™ Easy Preparation of ear notch samples 2.0). The presentation shows validation data for the ID Gene™ BVD / BD Triplex 2.0 and its companion kit for rapid direct lysis for ear notches.

## CONCLUSION

The ID Gene™ BVDV/BDV Triplex 2.0 kit offers:

- **Excellent performance:** high analytical sensitivity ; LDPCR = 5 copies / PCR and LDMMethod ≤ 25 TCID<sub>50</sub>/ ear notch when processed with the rapid direct lysis kit ID Gene™ Easy Preparation of Ear Notch Samples 2.0 (EZNOTCHV2).
- Test **pools of up to 25 samples**, thanks to high kit sensitivity
- **Optimal test reliability** thanks to endogenous and new best-in-class RNA-based exogenous internal controls which provide reliable results and efficacy, especially when associated with the direct lysis kit ID GENE™ Easy Preparation of ear notch samples 2.0
- **Compatible** with the direct lysis kit available in liquid format with economical **room temperature shipping** worldwide
- An included **positive amplification control (PAC)**. Sold separately : a positive extraction control calibrated to validate the entire PCR process

## ANALYSIS PROTOCOL

The ID Gene™ BVD/BD Triplex 2.0 (IDBVDV2) kit is a ready-to-use RT-qPCR kit assay detecting simultaneously BVD/BD virus, an endogenous and an exogenous internal positive controls.

| EXTRACTION              |  |
|-------------------------|--|
| Sample preparation      | Exogenous control and the Direct Lysis buffer were added to all individual ear notches samples           |
| Incubation              | 15 min at 100°C / Cooling step   |
| AMPLIFICATION           |  |
| Reaction condition      | 8µL of ready to use Mastermix was added to 5µL of each lysate from individual or up to 25 pooled samples |
| Amplification condition | Amplification were carried out using a rapid amplification program (65 min)                              |

| SAMPLE RESULT          | BVDV/BDV SIGNAL | ENDOGENOUS CONTROL SIGNAL | EXOGENOUS CONTROL SIGNAL | INTERPRETATION   |
|------------------------|-----------------|---------------------------|--------------------------|--|
| Positive               | +               | +/-                       | +/-                      | Valid; Animal detected as positive for BVDV/BDV<br>No signal for endogenous and/or exogenous control because of competition  |
| Negative               | -               | +                         | +                        | Valid; Animal detected as negative for BVDV/BDV  |
| Degraded sample        | -               | -                         | +                        | Not Valid; Analytical process successful, no signal for endogenous control due to poor sample quality – other sample analysis is required  |
| Partial PCR inhibition | -               | +                         | -                        | Not Valid; no signal for exogenous control or Cq exogenous on sample > Cq exogenous on negative extraction control + 3<br>Dilute the extracted DNA 10 times in nuclease-free water |

## RESULTS

### INCLUSIVITY

Analysis of 55 reference BVDV and BDV strains provided by the French Laboratory (LABOCEA Fougères) and by Friedrich Loeffler Institute (FLI, Gemany)

| BVDV/BD ISOLATE                | IDBVDV2 RESULTS | BVDV/BD ISOLATE                | IDBVDV2 RESULTS |
|--------------------------------|-----------------|--------------------------------|-----------------|
| FLI (Germany)                  | DETECTED        | 197-1 = BVD1 -b                | DETECTED        |
| LNK-BVD (France)               | DETECTED        | 207-1 = BVD1 -e                | DETECTED        |
| C strain FLI                   | DETECTED        | 9001 = BVD1 -e                 | DETECTED        |
| Eystrup91 FLI                  | DETECTED        | 09005-4 = BVD1 -b              | DETECTED        |
| Alfort187 FLI                  | DETECTED        | 09002-1 = BVD1 -e              | DETECTED        |
| Koslov1128 FLI                 | DETECTED        | 08185-1 = BVD1 -d              | DETECTED        |
| Brescia FLI                    | DETECTED        | 08158-5 = BVD1-a               | DETECTED        |
| Schweiz II FLI                 | DETECTED        | 08185-2 = BVD1 -l              | DETECTED        |
| Pader FLI                      | DETECTED        | 13035-4 = BVD 2                | DETECTED        |
| Bergen FLI                     | DETECTED        | 4544-01 v12                    | DETECTED        |
| D4886/82/RO FLI                | DETECTED        | 14011 = BVD cythopathogene     | DETECTED        |
| Uelzen FLI                     | DETECTED        | 07033-1 = BVD atypique PCR neg | DETECTED        |
| Spante FLI                     | DETECTED        | VdVB 14 - 12                   | DETECTED        |
| CongenitalTremor FLI           | DETECTED        | VdVB 10 - 5                    | DETECTED        |
| Kanagawa FLI                   | DETECTED        | VdVB 02/02/2007                | DETECTED        |
| Moredun FLI                    | DETECTED        |                                |                 |
| Rudolph FLI                    | DETECTED        |                                |                 |
| Gifhorn FLI                    | DETECTED        |                                |                 |
| Isard FLI                      | DETECTED        |                                |                 |
| NADL FLI                       | DETECTED        |                                |                 |
| Papltitz FLI                   | DETECTED        |                                |                 |
| PI809 FLI                      | DETECTED        |                                |                 |
| NC3807-1251/1 FLI              | DETECTED        |                                |                 |
| Egbert FLI                     | DETECTED        |                                |                 |
| BO806-17 FLI                   | DETECTED        |                                |                 |
| NC3807-8757 FLI                | DETECTED        |                                |                 |
| 8644 FLI                       | DETECTED        |                                |                 |
| Bure FLI                       | DETECTED        |                                |                 |
| Walter FLI                     | DETECTED        |                                |                 |
| PO1600 FLI                     | DETECTED        |                                |                 |
| Hobi FLI                       | DETECTED        |                                |                 |
| Giraffe FLI                    | DETECTED        |                                |                 |
| NCP-2508-FCS FLI               | DETECTED        |                                |                 |
| Bohni FLI                      | DETECTED        |                                |                 |
| 018901 = BVD1 e/cythopathogene | DETECTED        |                                |                 |
| 08207-3 = BVD1 -b              | DETECTED        |                                |                 |
| 163-5 = BVD1 -b                | DETECTED        |                                |                 |
| 08167-1 = BVD1 -b              | DETECTED        |                                |                 |
| 172-1 = BVD1 -e                | DETECTED        |                                |                 |
| 176-2 = BVD1 -b                | DETECTED        |                                |                 |

The ID Gene™ BVD/BD Triplex 2.0 detects all the reference strains tested, which include both BVDV (Genotype 1 and 2) and BDV viruses, as well as circulating viral isolates

### EXCLUSIVITY

The following 32 selected isolates were tested:

| SAMPLE NUMBER | STRAINS  | IDBVDV2 RESULTS |
|---------------|--|-----------------|
| 1-2           | Influenza (H1N1, H5N2)                               | Not detected    |
| 3-4           | Bovine Herpesvirus (1,4)                             | Not detected    |
| 5-8           | Bovine Pneumovirus (3,A,B)                           | Not detected    |
| 9-10          | PRRSV (EU and NA types)                              | Not detected    |
| 11-32         | 22 other viruses and bacteria from different species | Not detected    |


The ID Gene™ BVD/BD Triplex 2.0 kit did not show any cross-reactions with the 32 other pathogens tested

### OVERALL PERFORMANCE

| PCR CHARACTERISTICS  |  |
|--|--|
| Efficiency / R <sup>2</sup>  | 100.6 % / 0,98   |
| Limit of detection   | 5 copies/PCR with rapid amplification program                          |
| Robustness   | Unaffected by all parameters tested (temperature and volume of sample) |
| METHOD CHARACTERISTICS<br>(ID GENE™ EASY PREPARATION OF EAR NOTCH SAMPLES 2.0 : DIRECT LYSIS EXTRACTION) |  |
| Experimental limit of detection on ear notch sample spiked with BVDV strain (C24V Oregon)                | 25 TCID <sub>50</sub> /ear notch with rapid amplification program      |

### COMPARISON OF BOTH DIRECT LYSIS PROTOCOLS


The performance of the two direct lysis protocols were compared on 10 individual characterized field samples and 10 pools of 25 samples. The lysis protocol v1 (Direct Lysis v1) has a preliminary incubation step (70° C, 30min) compared to the new version of the lysis protocol (Direct Lysis v2). Direct Lysis v1 was combined with first version of IDBVD kit whereas the new Lysis protocol v2 was evaluated with IDBVDV2.



- The Direct Lysis v1 and Direct Lysis v2 show excellent agreement
- The Direct Lysis v2 gave lower Cq values than the Direct Lysis v1 (mean difference ΔCq = 3.54)
- One inhibiting positive pooled sample (Cq > 40) was not detected with Direct Lysis v1 but found positive with Direct Lysis v2 which ascertain its better ability to manage the inhibitors present in the lysates.

### COMPARISON WITH 2 PCR KITS

The performances of the new PCR kit were compared with the first version of ID GENE™ PCR kit (IDBVD) on 10 individual characterized field samples and 10 pooled to 25 samples. All individual or pooled samples were extracted with automated extraction method (with ID GENE™ MAG FAST Extraction Kit)



- New PCR kit (IDBVDV2) gave lower Cq values than the former version of PCR kit (IDBVD; mean difference ΔCq = 0.93)